

# Single High-pressure Six-way Switching Valve / Dual High-pressure Six-way Switching Valve

---



## Indispensable for Column Switching Analysis

Various analytical systems can be constructed by combining a great variety of components.

### Low Dead Volume Switching

The unique low-dead-volume valve head minimizes the influence of time lag and diffusion on separation in column switching analysis. (Dead volume: 0.8  $\mu$ L)

### Metal-free Structure

The flow channel is made of PEEK resin, making the valve effective for analyzing proteins of high adsorption to metals and physiologically active samples.

### Programming Function

A schedule can be programmed easily to automate complex column switching.

### Optimum Even for Complicated Multipurpose Channels

Either the single or dual type can be selected according to your needs.

### Integration to Other HPLC Systems

Even a start signal from another manufacturer's system can activate the self program, and so can be used as part of the system.



### Standalone or Networked Control

The valve allows programmed operation (standalone) by numerical keys without any special controller and has a function to communicate with Osaka Soda system controller program (PC).

### Example of Single Valve Application

Analysis of diazepam

Combination of single high-pressure six-way switching valve and CAPCELL PAK MF (column for direct introduction of biological samples) for automatic analysis of diazepam in human serum




### Primary Separation (A + B)

Column : CAPCELL PAK MF cartridge  
Mobile phase : 100 m mol/L Phosphate buffer/Acetonitrile = 98/2, pH7.0  
Flow rate : (A) 1.0 mL/min  
(B) 200  $\mu$ L/min  
Temperature : 35°C  
Injection volume : 80  $\mu$ L

### Secondary Separation (C)



Column : CAPCELL PAK C18 UG120  
2.0 mm i.d. x 250 mm  
Mobile phase : 100 m mol/L Phosphate buffer/  
Acetonitrile = 50/50, pH7.0  
Flow rate : 200  $\mu$ L/min  
Temperature : 35°C  
Detection : UV 230 nm



### Example of Dual Valve Application

Accurate gradient curve at low flow rate

Combination of dual high-pressure six-way switching valve and CAPCELL PAK MF (column for direct inspection of biological samples) for trace analysis of Chinese medicine and metabolites in mouse serum (Gradient system)




Source: Prof. Tashiro of the Showa Pharmaceutical University

Example of Dual Valve Application

## Rinsing system

The valve can also be built into a system for washing away sample carryover.



## Specifications

Product No.	Single high-pressure six-way switching valve: 3011 Dual high-pressure six-way switching valve: 3012
Product Name	Single high-pressure six-way switching valve/Dual high-pressure six-way switching valve
System	6-port, 2-position high-pressure switching
Number of Valves	3011: One 3012: Two
Max. Pressure	35 MPa
Channel Material	Polyetheretherketone (PEEK), ceramic
Control Mode	(1) Manual control (2) Control by self-program (0.00 to 999.99 min) (3) Control by external input (4) System control
Dimensions	120(W) x 230(H) x 491(D) mm
Weight	3011: About 8 kg 3012: About 10 kg