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

## Quinic acid

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
Quinic acid, a lower fatty acid, is contained in cranberry and coffee beans. It keeps urine weakly acidic, and prevents an urinary passage from infective disease. While it is highly hydrophilic and hardly retained in a reversed-phase column, quinic acid could be separated well with a NH<sub>2</sub> column based on weak anion-exchange chromatography. The following results were obtained after a simple pretreatment on coffee.

### 【HPLC Conditions 1】

1. キナ酸  
Quinic acid (M.W. 192.2)

### 【HPLC Conditions 2】

**【HPLC Conditions 1】**

Column : CAPCELL PAK NH<sub>2</sub> UG80 S5 ; 2.0 mm i.d. x 150 mm  
Mobile phase : 10 mmol/L Phosphate buffer (KH<sub>2</sub>PO<sub>4</sub> : Na<sub>2</sub>HPO<sub>4</sub> = 1 : 1 in molar ratio) / CH<sub>3</sub>OH = 98 / 2  
Flow rate : 200 μL/min  
Temperature : 40 °C  
Detection : PDA 210 nm  
Inj. vol. : 2 μL  
Sample dissolved in : Standard was dissolved in mobile phase.  
Ammonium hydroxide solution was added to coffee to make it at pH 9. A 0.5-mL portion of the treated coffee and 0.5 mL of chloroform were placed in a 1-mL vial. The resultant liquid was vigorously mixed and centrifuged. The chloroform layer was discarded. The addition and extraction with chloroform was repeated three times. The cleaned water layer was passed through a 0.45-μm filter and ten-fold diluted with the mobile phase.

**【HPLC Conditions 2】**

Column : CAPCELL PAK NH<sub>2</sub> UG80 S5 ; 2.0 mm i.d. x 250 mm  
(Two columns connected in series)  
Mobile phase : 10 mmol/L Phosphate buffer (KH<sub>2</sub>PO<sub>4</sub> : Na<sub>2</sub>HPO<sub>4</sub> = 1 : 1 in molar ratio) / CH<sub>3</sub>OH = 98 / 2  
Flow rate : 160 μL/min  
Temperature : 40 °C  
Detection : PDA 210 nm  
Inj. vol. : 2 μL  
Sample dissolved in : Sample prepared according to the same procedure as that in Conditions 1.  
※ 1 μg/mL = 1 ppm