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
Poeny belongs to Paeoniaceae. Its dried root is used as a Chinese medicine. The dried root contains 2% of paeoniflorin, a glycosilated monoterpene. Paeoniflorin is known to cause teleangiectasia and show anti-inflammatory effect.

The Japanese Pharmacopeia designates an HPLC method with a C18 column for its determination. A reference compound used for its system suitability test was switched from *p*-hydroxyacetophenone to albiflorin in the 15th Revision.

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
System suitability test

The Japanese Pharmacopeia 15th revision includes the following items in system suitability test of the compound.

1. Eluting order of albiflorin and, thereafter, paeoniflorin
2. Resolution between the above two compounds should be 2.5 or greater.
3. Paeoniflorin peak area of six replicates should show relative standard deviation of 1.5% or smaller.

The method with a CAPCELL PAK C₁₈ MGII S5 passed all of the above three items (1. passed, 2. resolution = 5.1, and 3. RSD = 0.32%).

Pretreatment : Peony root was cut into small pieces, and ground in a mortar with a pestle. 0.5 g of the ground peony root was dispersed in a mixed solvent (50 mL, water / methanol = 1 / 1). The dispersion was refluxed for 30 min. After cooling it down, the dispersion was filtered with filter paper. The residue was rinsed with the mixed solvent. The filtered liquid and the mixed solvent used for rinsing were added together and adjusted to 100 mL in a volumetric flask. A small amount of the solution was further filtered with a 0.2- μ m filter, and introduced to HPLC.