

**Headquarters**

12-18, Awaza 1-chome, Nishi-ku, Osaka 550-0011, Japan  
Telephone:+81-6-6110-1560 Facsimile:+81-6-6110-1603

January, 2020

Dear Valued Customers

**CAPCELL PAK MG Series Raw Material Change Schedule Information**

Osaka Soda Co., Ltd.  
Chemical Specialties Division  
Chromatography Sales Department

Thank you for using our high-quality CAPCELL PAK MG series HPLC columns.

As it was announced in December 2018, there was need to change the raw material for the CAPCELL PAK MG series columns.

Since we have secured 3 separate lots of the new material the change schedule is announced hereby. We apologize for any inconvenience caused by changing the raw materials.

**1. Shipping Date and Products**

- (1) Products : CAPCELL PAK MG/MG II /MG III /MG III-H
- (2) Change time : From July 1st, 2020

The method of identifying products before and after the change of raw materials will be described on the following page.

Should you have any questions or concerns, please inquiry to the following contact details.

Thank you for your understanding.

Osaka Soda Co., Ltd.  
Chemical Specialties Division  
Chromatography Sales Department  
12-18, Azawa 1-chome, Nishi-ku, Osaka 550-0011, Japan  
TEL: 06-6110-1598, FAX: 06-6110-1612  
E-mail: [silica@osaka-soda.co.jp](mailto:silica@osaka-soda.co.jp)

**Headquarters**

12-18, Awaza 1-chome, Nishi-ku, Osaka 550-0011, Japan  
 Telephone:+81-6-6110-1560 Facsimile:+81-6-6110-1603

## How to identify the raw material

### 1. Column number

Old raw material: ABCD10000 ~ ABCD49999

New raw material : ABCD50001 ~

### 2. Batch number

	S-3(3 $\mu$ m) Batch number		S-5(5 $\mu$ m) Batch number	
	Old	New	Old	New
MG	BMS- * *	BMO- * *	BS- * *	BO- * *
MG II	BMS II- * *	BMO II- * *	BS II- * *	BO II- * *
MG III / MG III-H	BMS III- * *	BMO III- * *	BS III- * *	BO III- * *

### 3. Label color

Old raw material (MG I/MG II)



New raw material



Old raw material (MG III/MG III-H)

