



1	Ketoprofen	2	Flurbiprofen	3	Indomethacin
	15 ppm		10 ppm		10 ppm
4	Ibuprofen				
	15 ppm				

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN / H ₂ O = 60 / 40, pH2.2 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 230 nm