



1	Ketoprofen	2	Flurbiprofen	3	Indomethacin
	15 µg/mL		10 µg/mL		10 µg/mL
4	Ibuprofen				
	15 µg/mL				

#### HPLC Conditions

Column	Supiorex ODS : 4.6 mm i.d.×250 mm
Mobile phase	CH <sub>3</sub> CN / H <sub>2</sub> O = 55 / 45, pH2.2 (H <sub>3</sub> PO <sub>4</sub> )
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 230 nm