



1	Primidone 15 ppm	2	Ethosuximide 250 ppm	3	Phenobarbital 10 ppm
4	Carbamazepine 25 ppm	5	Phenytoin 5 ppm		

#### HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	25mmol/L Phosphate buffer / CH <sub>3</sub> CN = 60 / 40
Flow rate	1.0 mL/min
Temperature	30 °C
Detect	UV 220 nm