



1	Primidone 15 µg/mL	2	Ethosuximide 250 µg/mL	3	Phenobarbital 10 µg/mL
4	Carbamazepine 2.5 µg/mL	5	Phenytoin 5 µg/mL		

HPLC Conditions

Column	Superiorex ODS : 4.6 mm i.d.×250 mm
Mobile phase	25mmol/L KH ₂ PO ₄ / CH ₃ CN / 25mmol/L Na ₂ HPO ₄ = 64 / 33 / 3, pH6.29
Flow rate	1.0 mL/min
Temperature	R.T.
Detect	UV 220 nm

