



- | | | |
|--|------------------|----------------------------|
| 1 Caffeine | 1 , cetaminophen | 2 Dihydrocodeine phosphate |
| 3 <i>dl</i> -Methylephedrine hydrochloride | 4 Noscapine | 5 Carbinoxamine maleate |
| 6 Bisibuthiamine | | |

HPLC Conditions

Column	Superiorex ODS : 4.6 mm i.d.×250 mm
Mobile phase	10mmol/L SDS, CH ₃ CN / H ₂ O = 60 / 40, pH2.2 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40 °C

Detect

UV 230 nm