



1	Scopolamine hydrobromide	0.2 mg/mL	2	Atropine sulfate	0.2 mg/mL
---	--------------------------	-----------	---	------------------	-----------

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	10mmol/L NaH ₂ PO ₄ / CH ₃ CN = 85 / 15, pH2.0
Flow rate	1.0 mL/min
Temperature	R.T.
Detect	UV 210 nm
Inj. vol.	10 μL