



1 Procaine hydrochloride

2 Caffeine(I.S.)

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG120 S5 : 4.6 mm i.d.×150 mm
Mobile phase	0.1vol% 1-Pentanesulfonic acid sodium salt, 50mmol/L KH ₂ PO ₄ (pH3.0, H ₃ PO ₄) / CH ₃ OH = 60 / 40
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254 nm
Inj. vol.	5 µL
Resolution	3.7