



- 1 Propyl *p*-hydroxybenzoate      2 Butyl *p*-hydroxybenzoate (I.S.)      3 Cortisone acetate

#### HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×150 mm
Mobile phase	H <sub>2</sub> O / CH <sub>3</sub> CN = 65 / 35
Flow rate	1.5 mL/min
Temperature	25 °C
Detect	UV 254 nm
Inj. vol.	10 µL
Resolution	12.4