



- 1 Ethyl *p*-hydroxybenzoate 2 Propyl *p*-hydroxybenzoate 3 Diclofenac sodium

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG120 S5 : 4.6 mm i.d.×150 mm
Mobile phase	CH ₃ OH / 0.12vol% CH ₃ COOH = 4 / 3
Flow rate	1.2 mL/min
Temperature	40 °C
Detect	UV 240 nm
Inj. vol.	10 µL
Resolution	5.4