



- | | | |
|-----------------------------------|---|------------------------------------|
| 1 Dexamethasone | 2 Propyl <i>p</i> -hydroxybenzoate (I.S.) | 3 Methyl <i>p</i> -hydroxybenzoate |
| 4 Ethyl <i>p</i> -hydroxybenzoate | | |

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	H ₂ O / CH ₃ CN = 2 / 1
Flow rate	1.8 mL/min
Temperature	25 °C
Detect	UV 254 nm
Inj. vol.	20 μL
Resolution	8.0