



1 Fluorometholone

2 Butyl *p*-hydroxybenzoate(I.S.)

HPLC Conditions

| | |
|--------------|---|
| Column | CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm |
| Mobile phase | CH ₃ OH / H ₂ O = 70 / 30 |
| Flow rate | 0.6 mL/min |
| Temperature | 35 °C |
| Detect | UV 254 nm |
| Inj. vol. | 20 µL |
| Resolution | 6.7 |