



1 Lidocaine hydrochloride

2 Benzophenone (I.S.)

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG120 S5 : 4.6 mm i.d.×150 mm
Mobile phase	2.88g SDS / 1000mL {20mmol/L phosphate buffer (pH3.0) / CH ₃ CN = 55 / 45}
Flow rate	2.3 mL/min
Temperature	25 °C
Detect	UV 254 nm
Inj. vol.	5 µL
Resolution	6.5