



1 Thymol (I.S.)

2 Propofol

HPLC Conditions

| | |
|--------------|--|
| Column | CAPCELL PAK C ₁₈ UG120 S3 : 4.6 mm i.d.×100 mm |
| Guard Column | CAPCELL PAK C ₁₈ UG120 S5 Cartridge : 4.0 mm i.d. × 10 mm |
| Mobile phase | CH ₃ CN / H ₂ O = 60 / 40 (v/v) |
| Flow rate | 800 μL/min |
| Temperature | 35 °C |
| Detect | FL Ex. 230 nm, Em. 340 nm |
| Inj. vol. | 20 μL |