



- | | | |
|-------------------|----------------|-----------------|
| 1 Ethoxybenzamide | 2 Phenytoin | 3 Ibuprofen |
| 4 Lidocaine HCl | 5 Progesterone | 6 Verapamil HCl |
| 7 Dibucaine HCl | | |

HPLC Conditions

Column	CAPCELL PAK C8 DD S3 : 4.6 mm i.d.×100 mm
Mobile phase	CH ₃ OH / 50mmol/L Phosphate buffer = 55 / 45
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 220 nm