



- | | | | | | |
|---|---------------|---|--------------|---|---------------|
| 1 | Ethenzamide | 2 | Phenytoin | 3 | Ibuprofen |
| 4 | Lidocaine | 5 | Progesterone | 6 | Verapamil HCl |
| 7 | Dibucaine HCl | | | | |

HPLC Conditions

Column	CAPCELL PAK C8 DD S5 : 4.6 mm i.d.×150 mm
Mobile phase	CH ₃ OH / 50mmol/L Phosphate buffer (pH=6.9) = 55 / 45
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 220 nm