



- | | | |
|--------------------|--------------------|-----------------|
| 1 Sulfadiazine | 2 Sulfathiazol | 3 Sulfamerazine |
| 4 Sulfadimidine | 5 Sulfamethoxazole | 6 Sulfisoxazole |
| 7 Sulfadimethoxine | | |

HPLC Conditions

Column	CAPCELL PAK C18 AQ S5 : 2.0 mm i.d.×20 mm
Mobile phase	(A) 0.1vol% HCOOH, H ₂ O (B) 0.1vol% HCOOH, CH ₃ OH B 0% - 50% (4.0min) - 0% (4.1min) Gradient
Flow rate	600 µL/min
Temperature	40 °C
Detect	UV 270 nm
Inj. vol.	2.0 µL