



1	Progesterone(I.S.)	22.5 µg/mL	2	Testosterone propionate	25 µg/mL
---	--------------------	------------	---	-------------------------	----------

HPLC Conditions

Column	CAPCELL PAK C18 MGII S5 : 4.6 mm i.d.×150 mm
Mobile phase	CH ₃ CN / H ₂ O = 7 / 3
Flow rate	1.0 mL/min
Temperature	35 °C
Detect	UV 241 nm
Inj. vol.	5 µL
Resolution	19.4