



- | | | |
|--------------------------|--------------------|------------------|
| 1 Sulfanilamide | 2 Sulfamethoxazole | 3 Sulfamethizole |
| 4 Sulfamonomethoxine | 5 Sulfamerazine | 6 Sulfapyridine |
| 7 Sulfamethoxypyridazine | 8 Sulfadimethoxine | 9 Sulfamethazine |
| 10 Sulfaquinoxaline | | |

HPLC Conditions

Column	CAPCELL PAK C8 DD S3 : 4.6 mm i.d.×250 mm
Mobile phase	20mmol/L Phosphate buffer / CH ₃ OH = 83 / 17
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 270 nm