



1	Sulfamerazine	2	Sulfadimethoxine	3	Sulfaquinoxaline
	10 µg/mL		10 µg/mL		10 µg/mL

#### HPLC Conditions

Column	CAPCELL PAK MF Ph-1 : 4.6 mm i.d.×150 mm
Mobile phase	100mmol/L Phosphate buffer / CH <sub>3</sub> CN = 90 / 10, pH6.98
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254 nm
Inj. vol.	20 µL