



- | | | | | | |
|---|----------------------|---|----------------------|---|----------------------|
| 1 | Ibuprofen metabolism | 2 | Ibuprofen metabolism | 3 | Ibuprofen metabolism |
| 4 | Ibuprofen metabolism | 5 | Ibuprofen metabolism | 6 | Ibuprofen |

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 1.5 mm i.d.×250 mm
Mobile phase	50mmol/L Potassium dihydrogen phosphate pH2.45 / CH ₃ CN Gradient
Flow rate	100 µL/min
Temperature	40 °C
Detect	UV 210 nm