



- |   |                          |   |             |   |                     |
|---|--------------------------|---|-------------|---|---------------------|
| 1 | Epigallocatechin         | 2 | Catechin    | 3 | Caffeine            |
| 4 | Epigallocatechin gallate | 5 | Epicatechin | 6 | Epicatechin gallate |
| 7 | Catechin gallate         |   |             |   |                     |

#### HPLC Conditions

Column	CAPCELL PAK C18 UG120 S3 : 4.6 mm i.d.×100 mm
Mobile phase	0.5vol% H <sub>3</sub> PO <sub>4</sub> , H <sub>2</sub> O / CH <sub>3</sub> OH = 82 / 18
Flow rate	0.8 mL/min
Temperature	40 °C
Detect	UV 280 nm