



1	Benzoic acid	2	Saccharin	3	Sorbic acid
	100 ppm		100 ppm		100 ppm

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	{25mmol/L KH ₂ PO ₄ , 25mmol/L Na ₂ HPO ₄ } / CH ₃ OH = 90 / 10
Flow rate	0.7 mL/min
Temperature	50 °C
Detect	UV 225 nm
Inj. vol.	5 µL