



1	Rutin	2	Myricetin	3	Quercetin
	200 ppm		200 ppm		200 ppm
4	Geraldol				
	200 ppm				

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG120 S3 : 1.5 mm i.d.×150 mm
Mobile phase	H ₂ O / CH ₃ CN / 0.5mol/L Oxalic acid = 80 / 20 / 2
Flow rate	200 µL/min
Temperature	40 °C
Detect	UV 350 nm
Inj. vol.	1 µL