



- | | | |
|------------|-------------------|--------------|
| 1 Glycerin | 2 Pentaerythritol | 3 Erythritol |
| 4 Xylitol | 5 Sorbitol | 6 Mannitol |

HPLC Conditions

Column	CAPCELL PAK NH ₂ UG80 S5 : 1.5 mm i.d.×250 mm
Mobile phase	CH ₃ CN / H ₂ O = 85 / 15
Flow rate	100 µL/min
Temperature	40 °C
Detect	RI
Inj. vol.	1 µL