



1	Lauric acid (C12: 0) 1 mg/mL	2	Linolenic acid (C18: 3) 0.1 mg/mL	3	Myristic acid (C14: 0) 1 mg/mL
4	Linoleic acid (C18: 2) 0.4 mg/mL	5	Palmitic acid (C16: 0) 3 mg/mL	6	Oleic acid (C18: 1) 1 mg/mL
7	Stearic acid (C18: 0) 5 mg/mL				

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S3 : 4.6 mm i.d.×100 mm
Mobile phase	CH ₃ CN / 0.1vol% H ₃ PO ₄ = 80 / 20
Flow rate	1.0 mL/min
Temperature	R.T.
Detect	UV 210 nm
Inj. vol.	20 µL