



- |   |                            |   |                 |   |                 |
|---|----------------------------|---|-----------------|---|-----------------|
| 1 | Ethyl $\gamma$ -linolenate | 2 | Ethyl myristate | 3 | Ethyl linoleate |
| 4 | Ethyl palmitate<br>oleate  | 5 | Ethyl stearate  |   |                 |

#### HPLC Conditions 1

Column            Superiorex ODS : 4.6 mm i.d.×250 mm  
 Mobile phase    CH<sub>3</sub>OH  
 Flow rate        1.0 mL/min  
 Detect            RI  
 Inj. vol.         10  $\mu$ L

#### HPLC Conditions 2

Column            Superiorex ODS : 20 mm i.d.×250 mm  
 Mobile phase    CH<sub>3</sub>OH  
 Flow rate        10 mL/min  
 Detect            RI  
 Inj. vol.         100  $\mu$ L (83.5 mg)