



1 Sorbic acid

#### HPLC Conditions

Column	CAPCELL PAK C <sub>18</sub> MG S5 : 1.5 mm i.d.×250 mm
Mobile phase	(20mmol/L Tetrabutylammonium chloride, 20mmol/L KH <sub>2</sub> PO <sub>4</sub> , 20mmol/L Na <sub>2</sub> HPO <sub>4</sub> CH <sub>3</sub> CN = 70 / 30
Flow rate	100 µL/min
Temperature	40 °C
Detect	UV 265 nm
Inj. vol.	1 µL