



Benzalkonium chloride
1
Standard 1 mg/mL each

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN / {0.1mmol/L NH ₄ H ₂ PO ₄ , 0.1vol% H ₃ PO ₄ , 0.1vol% Triethylamine, pH2.5} = 70 / 30
Flow rate	1.0 mL/min
Temperature	50 °C
Detect	UV 254 nm
Inj. vol.	8 µL