



1 Benzalkonium chloride

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 1.5 mm i.d.×250 mm
Mobile phase	CH ₃ CN / {0.1mmol/L NH ₄ H ₂ PO ₄ , 0.1vol% H ₃ PO ₄ , 0.1vol% Triethylamine, pH2.5} = 70 / 30
Flow rate	110 μL/min
Temperature	50 °C
Detect	UV 254 nm
Inj. vol.	0.5 μL