



- | | | |
|-------------------------|---|-----------|
| 1 Thioglycolic acid | 2 Cysteine | 3 Cystine |
| 4 Dithiodiglycolic acid | 5 <i>p</i> -Hydroxybenzoic acid @
(I.S.) | |

HPLC Conditions

Column	CAPCELL PAK C18 AG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	3.5mmol/L Sodium heptanesulfonate, 0.1vol% H ₃ PO ₄ / CH ₃ CN = 95 / 5
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 210 nm
Inj. vol.	10 μL