



Stearyl Glycyrrhetinate

1

Standard 25 ppm

Product 20000 ppm

HPLC Conditions 1

Column CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm

Mobile phase CH₃OH

Flow rate 1.0 mL/min

Temperature R.T.

Detect UV 254 nm

Inj. vol. 20 μL

HPLC Conditions 2

Column CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm

Mobile phase	CH ₃ CN
Flow rate	1.0 mL/min
Temperature	45 °C
Detect	UV 254 nm
Inj. vol.	20 μL