



1	Isopropyl-paraben 100 ppm	2	<i>n</i> -Propylparaben 100 ppm	3	Isobutyl-paraben 100 ppm
4	<i>n</i> -Butylparaben 100 ppm				

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG80 S5 : 4.6 mm i.d.×150 mm
Mobile phase	CH ₃ OH / H ₂ O = 60 / 40
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254 nm
Inj. vol.	2 μL