



1	Resorcinol 1000 ppm	2	Swertiamarin 200 ppm	3	<i>p</i> -Hydroxybenzoic acid 100 ppm
---	------------------------	---	-------------------------	---	--

HPLC Conditions

Column	CAPCELL PAK C18 MG S5 : 4.6 mm i.d.×250 mm
Mobile phase	H ₂ O / CH ₃ CN = 87 / 13 (v/v), pH2.5 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 238 nm
Inj. vol.	5 μL