



- | | | |
|---|---------------------------------|---------------------|
| 1 Biotin (VH) | 2 Riboflavin (VB ₂) | 3 Nicotinic acid |
| 4 Pyridoxine hydrochloride (VB ₆) | 5 Folic acid (VBc) | 6 Nicotinamide (PP) |

HPLC Conditions

Column	CAPCELL PAK NH ₂ UG80 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN / 25mmol/L KH ₂ PO ₄ = 80 / 20, pH2.0 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 220 nm
Inj. vol.	20 µL