



1	Biotin (VH) 80 µg/mL	2	Riboflavin (VB ₂) 80 µg/mL	3	Ascorbic acid (VC) 100 µg/mL
4	Nicotinic acid 20 µg/mL	5	Folic acid (VBc) 160 µg/mL	6	Nicotinamide (PP) 20 µg/mL

HPLC Conditions

Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN / 25mmol/L KH ₂ PO ₄ = 80 / 20, pH2.0 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 210 nm
Inj. vol.	10 µL

