



1	α -Tocopherol	2	β -Tocopherol	3	γ -Tocopherol
	125 ppm		125 ppm		125 ppm
4	δ -Tocopherol				
	125 ppm				

HPLC Conditions

Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d.×250 mm
Mobile phase	<i>n</i> -Hexane / Ethyl acetate = 90 / 10
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 295 nm
Inj. vol.	5 μ L

