



- | | | |
|--------------------------|-------------|--------------------------|
| 1 Vitamin K ₃ | 2 Vitamin A | 3 Vitamin D ₂ |
| 4 Vitamin D ₃ | 5 Vitamin E | 6 Vitamin E acetate |
| 7 Vitamin K ₁ | | |

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN / CH ₃ OH = 80 / 20
Flow rate	1.0 mL/min
Temperature	35 °C
Detect	UV 254 nm