



- |   |            |   |                   |   |            |
|---|------------|---|-------------------|---|------------|
| 1 | Vitamin K3 | 2 | Vitamin A acetate | 3 | Vitamin D2 |
| 4 | Vitamin D3 | 5 | Vitamin E acetate | 6 | Vitamin K1 |

#### HPLC Conditions

Column	CAPCELL PAK C18 AG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH <sub>3</sub> CN / CH <sub>3</sub> OH = 80 / 20
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254 nm