



1	Vitamin K3 2 µg/mL	2	Vitamin A acetate 40 µg/mL	3	Vitamin D2 22 µg/mL
4	Vitamin D3 24 µg/mL	5	Vitamin E acetate 130 µg/mL	6	Vitamin K1 40 µg/mL

HPLC Conditions

Column	Supiorex ODS : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN
Flow rate	1.5 mL/min
Temperature	40 °C
Detect	UV 254 nm
Inj. vol.	20 µL