



1	Pyridoxal 5'-phosphate	2	Pyridoxal	3	Pyridoxamine
	50 ppm		50 ppm		50 ppm
4	Pyridoxine				
	50 ppm				

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 1.5 mm i.d.×150 mm
Mobile phase	{50mmol/L KH ₂ PO ₄ , 5mmol/L Sodium octansulfonate} / CH ₃ OH = 80 / 20
Flow rate	100 µL/min
Temperature	30 °C
Detect	UV 230 nm
Inj. vol.	1 µL