



1	Vitamin K3 30 µg/mL	2	Vitamin K2 150 µg/mL	3	Vitamin E 750 µg/mL
4	Vitamin E acetate 30 µg/mL	5	Vitamin K1 30 µg/mL		

HPLC Conditions

Column	CAPCELL PAK C18 AQ S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN / CH ₃ OH = 50 / 50
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254/284/254 nm
Inj. vol.	10 µL