



1	Vitamin K3 7.5 µg/mL	2	Vitamin K2 50 µg/mL	3	Vitamin D2 75 µg/mL
4	Vitamin D3 50 µg/mL	5	Vitamin E 500 µg/mL	6	Vitamin E acetate 800 µg/mL
7	Vitamin K1 100 µg/mL				

HPLC Conditions

Column	CAPCELL PAK C18 ACR S3 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ CN / CH ₃ OH = 80 / 20
Flow rate	1.0 mL/min
Temperature	35 °C
Detect	UV 254 nm, 284 nm
Inj. vol.	10 µL