



1	Glucose	2	Maltose	3	Maltotriose
	1000 µg/mL		1000 µg/mL		1000 µg/mL
4	Maltotetraose	5	Maltopentaose	6	Maltohexaose
	1000 µg/mL		1000 µg/mL		1000 µg/mL
7	Maltoheptaose				
	1000 µg/mL				

HPLC Conditions

Column	CAPCELL PAK C18 AQ S5 : 4.6 mm i.d.×250 mm
Mobile phase	H ₂ O
Flow rate	1.0 mL/min
Temperature	35 °C
Detect	RI
Inj. vol.	20 µL