



1	Chondro δ Di-0S	2	Chondro δ Di-6S	3	Chondro δ Di-4S
	0.6 μ g/mL		0.6 μ g/mL		0.6 μ g/mL

HPLC Conditions

Column	CAPCELL PAK NH ₂ UG80 S5 : 1.5 mm i.d.×250 mm
Mobile phase	0.25mol/L KH ₂ PO ₄
Flow rate	100 μ L/min
Temperature	40 °C
Detect	UV 232 nm
Inj. vol.	10 μ L