



- | | | |
|-------------------|-------------------|------------------|
| 1 Saikosaponin c | 2 Saikosaponin f | 3 Saikosaponin a |
| 4 Saikosaponin b2 | 5 Saikosaponin b1 | 6 Saikosaponin d |

HPLC Conditions 1

Column CAPCELL PAK C₁₈ AG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase H₂O / CH₃CN = 50 / 50
Flow rate 1.0 mL/min
Temperature 40 °C
Detect UV 210 nm

HPLC Conditions 2

Column CAPCELL PAK C₁₈ SG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase H₂O / CH₃CN / CH₃OH = 45 / 40 / 15
Flow rate 1.0 mL/min
Temperature 40 °C
Detect UV 210 nm
Inj. vol. 5 µL