



1 Saikosaponin a

2 Saikosaponin d

HPLC Conditions 1

Column CAPCELL PAK C18 UG80 S5 : 4.6 mm i.d.×250 mm

Mobile phase CH₃OH / H₂O = 67 / 33

Flow rate 1.0 mL/min

Temperature 40 °C

Detect UV 210 nm

Inj. vol. 10 µL

HPLC Conditions 2

Column CAPCELL PAK C18 UG80 S20 : 4.6 mm i.d.×250 mm

Mobile phase CH₃OH / H₂O = 67 / 33

Flow rate 1.0 mL/min

Temperature 40 °C

Detect UV 210 nm

Inj. vol. 10 µL